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Ackreditering

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Referenser

1. Kasper, C. K. Laboratory tests for factor VIII inhibitors, their variation, significance and interpretation. Blood Coagulation and Fibrinolysis. 1991; 2(suppl. 1): 7-10).
2. Verbruggen, B. et al. The Nijmegen modification of the Bethesda assay for factor VIII:C

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MA

-
- inhibitors: Improved specificity and reliability. *Thromb Haemost* 1995; 73: 241-251.
3. Giles, A. et al. A detailed comparison of the performance of the Standard versus the Nijmegen modification of the Bethesda assay in detecting factor VIII:C inhibitors in the haemophilia A population of Canada. *Thromb Haemost* 1998; 79: 872-875.
4. Bipacksedel till Coatest[®] Factor VIII, Chromogenix, Mölndal.
5. Meijer, Piet. Standardization of the Factor VIII Inhibitor Assay. ECAT Foundation, May 2015.
6. Instrumenthandhavande Atellica Coag (Siemens Healthcare Diagnostics)
7. Miller, C.H. Laboratory testing for factor VIII and IX inhibitors in haemophilia: A review. *Haemophilia*. 2018 March; 24(2): 186-197.
8. Verbruggen, B. et al. Detecting and quantifying acquired functional inhibitors in hemostasis. Chapter 12 in *Quality in Laboratory Hemostasis and Thrombosis*, Second edition.