P-Triglycerider på Atellica (NPU04094)

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P-TG på Atellica (NPU04094)

Referenser

1. Nilsson-Ehle P, red. Laurells Klinisk kemi i praktisk medicin. Lund: Studentlitteratur 2003, 8:e upplagan sid 327-56.
2. Nordestgaard B G et al. Fasting Is Not Routinely Required for Determination of a Lipid Profile: Clinical and Laboratory Implications Including Flagging at Desirable Concentration Cutpoints- A Joint Consensus Statement from the European Atherosclerosis Society and European Federation of Clinical Chemistry and Laboratory Medicine. Clin Chem 2016;62:930-46.
3. Att förebygga aterosklerotisk hjärt-kärlsjukdom med läkemedel - behandlingsrekommendation. Information från Läkemedelsverket 5:2014.
4. 2019 ESC/EAS Guidelines for the management of dyslipidaemias: lipid modification to reduce cardiovascular risk. The Task Force for the management of dyslipidaemias of the European Society of Cardiology (ESC) and European Atherosclerosis Society (EAS). Eur Heart J 2020;1:111-88.
5. Colantonio DA et al. Closing the gaps in pediatric laboratory reference intervals: a CALIPER database of 40 biochemical markers in a healthy and multiethnic population of children. Clin Chem 2012;58:854-68.
6. Ridefelt P et al. Pediatric reference intervals for general clinical chemistry components - merging of studies from Denmark and Sweden. SJCLI 2018;78:365-372.
7. Rustad P et al. The Nordic Reference Interval Project 2000: recommended reference intervals for 25 common biochemical properties. SJCLI 2004;64:271-284.
8. Siemens produktblad: Atellica CH Triglycerides_2 (Trig_2) Rev. 04, 2022-03.
9. Instrumenthandhavande Atellica 20-79.
10. Atellica analysdata 20-139.
11. ABC Analyshantering 20-65.