Pt-Svett-test (NPU17183)

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Referenser

1. Davis PB. Cystic fibrosis since 1938. *Am J Crit Care Med.* 2006; 173:475-482.
2. Zirbes J, Milla CE. Cystic fibrosis related diabetes. *Paediatr Respir Rev* 2009; 10(3):118-123.
3. Rommens JM et al. Identification of the cystic fibrosis gene: Chromosome walking and jumping. *Science* 1989; 245: 1059-65.
4. Riordan JR et al. Identification of the cystic fibrosis gene:cloning and characterization of complementary DNA. *Science* 1989; 245: 1066-73.
5. Kerem BS et al. Identification of the cystic fibrosis gene: genetic analysis. *Science* 1989; 245: 1073-80.
6. Paranjape SM, Zeitlin PL. Atypical cystic fibrosis and CFTR-related diseases. *Clinic Rev Allergy* 2008; 35:116-123.
7. Sant'Agnes di, PA, Darling RC, Perera GA, Shea E. Abnormal electrolyte composition of sweat i cystic fibrosis of the pancreas. *Pediatrics* 1953; 12:549-563.
8. Gibson L E, Cooke RE. A test for concentration of electrolytes in sweat in cystic fibrosis of the pancreas, utilising pilocarpine by iontophoresis. *Pediatrics* 1959; 23: 545-549.
9. Farrell PM, et al.Diagnosis of Cystic Fibrosis: Consensus Guidelines from the Cystic Fibrosis Foundation. *J Pediatrics* 2017; 181S:S4-S15.
10. CLSI. Sweat Testing: Sample Collection and Quantitative Chloride Analysis; approved guideline. 4th edition. CLSI guideline C34. Wayne PA: Clinical and Laboratory Standards Institute; 2019.
11. Cole DEC, Boucher MJ. Use of a new sample-collection device (Macroduct™) in anion analysis of human sweat. *Clin Chem* 1986; 32(7): 1375-1378.
12. Macroduct Advanced Sweat Collection System Model 3710 SYS. User´s Manual. 2021.
13. Dokument från Arbetsgruppen för Cystisk fibros upplaga [190707](#). Svetttest. Rekommendationer avseende utförande och tolkning. 2019.
14. Kloosterboer M et al. Clarification of laboratory and clinical variables that influence cystic fibrosis newborn screening with initial analysis of immunoreactive trypsinogen. *Pediatrics* 2009; 123:e338-e346
15. Sharp JK, Rock MJ. Newborn screening for cystic fibrosis. *Clinic Rev Allerg Immunol* 2008; 35:107-115.
16. Sands D et al. Bilateral sweat tests with two different methods as a part of cystic fibrosis newborn screening (CF NBS) protocol and additional quality control. *Folia Histochemica cytobiologica* 2010; 48(3):358-365.
17. Hammond KB et al. Clinical evaluation of the macroduct collection system and conductivity analyzer in the diagnosis of cystic fibrosis. *J Pediatr* 1994; 124(2):255-260.
18. Koerbin G et al. Total intra-individual variation in sweat sodium and chloride concentrations for the diagnosis of cystic fibrosis. *Clin Chim Acta* 2008; 393:128-129.
19. [Cystisk fibros - Socialstyrelsen \(2023-06-15\)](#)